# 3M<sup>™</sup> Products and Services for Medical Sterilization Processes





### 3M Medical Sterilization Specialties

3M Medical Sterilization Specialties is a global leader in the manufacture and supply of biological and chemical indicators for monitoring sterilization processes. 3M has led the way in developing innovative technologies and products including the first self-contained biological indicator and the unique 3M<sup>™</sup> Attest<sup>™</sup> Rapid Readout Biological Indicators for Ethylene Oxide and Steam sterilization processes.

In an increasingly competitive global manufacturing environment, 3M Attest Rapid Readout Biological Indicators provide the fastest (3 or 4 hours) biological indicator result available today providing a competitive advantage by maximizing inventory turns and process throughput leading to quantifiable economic advantages in manufacturing efficiency.

#### Contents

- 06 Monitoring Products for Ethylene Oxide (EtO) Sterilization
- 11 Sterilization Equipment
- 14 Monitoring Products for Steam Sterilization
- 20 Monitoring Products for Vapor Hydrogen Peroxide Sterilization
- 21 Monitoring Products for Dry Heat Sterilization
- 23 Other Sterilization Monitoring Products
- 25 Performance Standards References
- 26 Index

### 3M Health Care Profile and History

3M is a diversified technology company known worldwide for innovative and reliable products. With annual sales of more than \$4 billion, 3M Health Care is one of the company's largest business. The 3M Health Care family of businesses serves a varied marketplace including medical, dental, pharmaceutical, animal health, and health information markets. The extended 3M Health Care family also serves the personal care and food safety markets. We apply our broad experience of working with blue chip medical manufacturing companies to deliver our range of innovative sterilization assurance monitoring products to the healthcare industry.

#### 3M Health Care Philosophy

3M Health Care is committed to providing the health care manufacturing community with innovative solutions, based on 3M diversified technologies, that help improve your business efficiency.

### 3M Medical Sterilization Specialties

3M Medical Sterilization Specialties is a business within 3M Health Care. 3M has provided sterilization assurance monitoring products and hardware to the healthcare industry for over 25 years. We recognize that many of our sterilization assurance products have an application in industrial sterilization. Our business is focused specifically on the needs of Industrial Sterilization customers. As a medical device and pharmaceutical manufacturer ourselves, we have the same rigorous compliance and productivity requirements that drive your business.

3M Medical Sterilization Specialties has the unparalleled benefit of the wider 3M organization to back our products. Access to state-of-the-art R&D laboratories and Corporate Technology Centers combined with our renowned capabilities in technical excellence and 3M Innovation, provide a unique blend of skills to bring to the market highly innovative solutions for the health care manufacturing environment.

3M has over 25 separate corporate technology platforms to draw upon to support our customer's product needs. 3M Medical Sterilization Specialties products incorporate technologies from multiple 3M technology platforms including:

Adhesives
Biological sciences
Microbiology
Electromechanical systems
Fluorochemicals
Reactive Inks & pigments
Specialty films and materials

It is our access to these technology platforms and 3M Corporate technical expertise that provides 3M Medical Sterilization Specialties with the ability to provide unmatched product solutions.



### Technical Service and Product Support

Our Technical Service group is an integral part of our laboratory organization. Staffed with dedicated, experienced personnel with backgrounds in microbiology, biochemistry, food services, medical technology and other related disciplines, the group's primary responsibility is to provide direct support to 3M customers by responding to product and process related issues in a timely manner. The Technical Service team not only addresses product related questions, but can also help provide support to your cycle development and process validation projects, drawing from our extensive experience in sterilization.

Our Technical Service team is recognized as experts in their chosen field and keep current by participating within industry standards setting organization such as ISO, AAMI and other national standards setting groups. Participation in specific congresses around the world is an integral part of the team's activities.

3M Health Care is a global company with subsidiary organizations around the world. Our international operations are based in St. Paul, Minnesota where our primary Technical Service resources are located. However in recognition of the fact that our customers are also global companies, we have created a support structure in each of our regional organizations with technical expertise that can address issues in your local language at a convenient time. Our regional organizations have instant access to our product experts in St. Paul ensuring a rapid response to customer support requests.



3M Tech Line (US and Canada only) 1 800 441 1922 option 2

### **Biological Indicator Systems**



#### 3M™ Attest Rapid Readout Biological Indicator for Ethylene Oxide

3M Attest Rapid Readout biological indicator system for use with ethylene oxide (EtO) sterilization processing is a proprietary technology that provides the fastest biological indicator result available. 3M Attest Rapid Readout provides results in as little as four hours post-processing.

- Highly resistant, standardized population of Bacillus atrophaeus spores for EtO
- Provides all of the benefits of self-contained technology plus results in four hours
- · Reduced incubation can improve inventory levels and increase turns
- · Meets FDA classification as a biological indicator
- Complies with ISO and ANSI/AAMI/ISO 11138:2006, Parts 1 and 2

Cat. No.	Description	Packaging	Shipping
3M™ Attest	™ Rapid Readout Biological Indicator for Ethylen	e Oxide	
1294-S	Biological Indicator for Ethylene Oxide	300 each/box	2 boxes/case



#### 3M™ Attest™ 293G and 290G Auto-readers

3M 293G and 290G Auto-readers are designed to incubate and automatically read the fluorescence of the Attest 1294-S Rapid Readout Biological Indicator in 4 hours. The large capacity 3M Attest 293G Auto-reader is designed specifically for industrial sterilization processing.

- Automatic readout with visual display of indicator result
- The 293G has a 36 well capacity and the 290G has a 12 well capacity
- Developed with significant input from leading medical device manufacturers
- Can help streamline your sterilization monitoring process
- Incubates at  $37^{\circ}\text{C} \pm 2^{\circ}\text{C}$



Cat. No.	Description	Packaging	Shipping
3M™ Attest	t™ Auto-readers		
293G	Auto-reader for Ethylene Oxide Large capacity 3M Attest Auto-reader For use with 3M Attest 1294-S Rapid Readout for Ethylene Oxide Biological Indicator 36 well capacity	1 unit/box	1 box/case
290G	Auto-reader for Ethylene Oxide For use with 3M Attest 1294-S Rapid Readout for Ethylene Oxide Biological Indicator 12 well capacity	1 unit/box	1 box/case

### Conventional, Self-contained Biological Indicators

#### 3M™ Attest™ Biological Indicators for Ethylene Oxide (EtO) Sterilization

3M developed the first self-contained biological indicators to reduce contamination from improper transfer technique.

- Highly resistant, standardized population of Bacillus atrophaeus spores for EtO
- Conventional self-contained technology benefits:
  - Simplistic culturing
  - Minimizes false positives due to contamination
  - Cost-effective
- Routinely documented for incubation times of less than seven days (typically 48 hours)
- Complies with ISO and ANSI/AAMI/ISO 11138:2006, Parts 1 and 2

Cat. No.	Description	Packaging	Shipping
3M™ Attest	™ Biological Indicators for Ethylene Oxide (EtO	) Sterilization	
1264-S	EtO Biological Indicator (Green cap)	300 each/box	2 boxes/case

#### 3M™ Attest™ Biological Indicator Incubator

Attest Incubators provide optimal conditions for accurate readout of Attest 1264-S Biological Indicators.

- Convenient, 28-vial capacity visible from above
- Dry-block design pre-set to proper incubation temperature

Cat. No.	Description	Packaging	Shipping
3M™ Attes	t <sup>™</sup> Biological Indicator Incubator for 1264-S		
127	28-vial (2-tier), 37°C (Et0), 120 nominal voltage required, 90–132V acceptable (U.S.)	1 each/box	1 box/case
129	28-vial (2-tier), 37°C (Et0) 220/240 volt	1 each/box	1 box/case





### Spore Strips and Media



### 3M<sup>™</sup> Attest<sup>™</sup> Spore Strips and Culture Media for Ethylene Oxide or Dry Heat

Spore strips for industrial sterilization applications

- Highly resistant, standardized population of Bacillus atrophaeus spores
- Packaged in peel-open glassine paper pouch
- Comply with ISO and ANSI/AAMI/ISO 11138:2006, Parts 1 and 2 (EtO)

#### **Culture Media**

Tryptic Soy Broth (TSB) for use with 3M Attest Spore Strips for EtO or Dry Heat

- Meets USP recommendation for soybean casein digest broth
- Meets NCCLS sterility and growth promotion guidelines

Cat. No.	Description	Packaging	Shipping
3M™ Attest	Spore Strips for EtO or Dry Heat		
1263-S	Biological indicator spore test strip for EtO or Dry Heat	100 each/box	4 boxes/case
Culture Med	lia Test Kit		
1263-M	25 1263-S spore strips and 25 TSB media tubes/box	25 of each/box	4 boxes/case

### **Chemical Indicators**

#### **Class 5 Chemical Integrating Indicators for Ethylene Oxide**

For use in parametric release processes when independent monitoring information is indicated.

- · Indicates sterilization critical parameters have been met
- 1244 Moving front technology provides easy-to-read accept/reject results
- For use with EtO/HCFC gas mixtures (1244 A)
- 1251 color change is red to green when exposed to ethylene oxide sterilization conditions
- Comply with ISO and ANSI/AAMI/ISO 11140:2005, Part 1

Cat. No.	Description	Packaging	Shipping
3M™ Comp	ly™ (SteriGage™) Ethylene Oxide (Et0) Chemical	Integrating Indicators	
1244A	EtO Chemical Integrator, moving front, for gas mixture sterilizers, 5,1 cm x 1,9 cm	500 each/bag ı (2" x 3/4")	2 bags/case
3M™ Comp	ly™ Ethylene Oxide (EtO) Chemical Indicator Stri	p	
1251	EtO Chemical Indicator Strip	240 each/box	1 box/case





#### Process Indicators and Accessories



#### **Indicator Tapes for EtO Sterilization**

A reliable ethylene oxide gas indicator tape that provides immediate identification of processed goods. Chemical indicator lines will turn red when exposed to an EtO sterilization process.

· Reliable adhesive adheres to packaging material

Cat. No.	Description	Packaging	Shipping
2M™ Comp		limation#	
SIVI COITIP	ly™ Indicator Tape for Ethylene Oxide (Et0) Steri	IIZation*	

<sup>\*</sup>Contains Dry Natural Rubber





#### **Indicator Tape Dispensers**

**Definite length** 

M920

Cat. No.	Description	Packaging	Shipping
C22 Heavy	Duty Tape Dispenser		
two rolls of	eavy duty tape dispenser is made of rugg f 3/4-inch (1 cm) wide indicator tape. ted base	ed, high-impact, durable plas	stic. Capable of handling
C22	Heavy duty, holds 2 rolls	_	2 each/case
M920 Defir	nite-length Tape Dispenser		
	definite-length tape dispenser adjusts to	dispense tape at a set lengtl	n.
	ngth up to 4 inches (10 cm) one indicator tape roll up to 1 inch (2.5 c	m) in width	

6 each/case

### Ethylene Oxide Sterilizers/Aerators

#### **Ethylene Oxide Sterilizers/Aerators**

3M Steri-Vac sterilizers provide dual cycle (37°C and 55°C), 100% EO gas sterilization for small-chamber, inhouse industrial sterilization.

- Supplied with user adjustable process parameters
- Single dose 100% EtO cartridge
- · Engineered for safety
- Built-in aeration cycle
- Effective low temperature sterilization process for industrial medical, pharmaceutical and veterinary applications



#### 3M™ Steri-Vac™ 8XL Sterilizer/Aerator

In addition to all of the features of the 5XL, the 8XL provides a 200 L (7.9 cu. ft.) chamber. The relative humidity feature continuously monitors, maintains and documents the proper humidity in the sterilization chamber during the precondition phase.

8XL Sterilizer/aerator 1 unit/case

Exterior: 179 cm (H) x 94 cm (W) x 109 cm (D)

70" (H) x 37" (W) x 43" (D) Chamber: 46 cm (H) x 51 cm (W) x 97 cm (D) 18-1/8" (H) x 20" (W) x 38" (D)

#### 3M™ Steri-Vac™ 5XL Sterilizer/Aerator

The 5XL sterilizer features a 130 L (4.8 cu. ft.) chamber. The validated computer controlled sensors ensure accuracy by automatically abating the cycle and displaying an error code if a malfunction is detected. Relative humidity detection and documentation is available.

5XL Sterilizer/aerator

Exterior: 59 cm (H) x 76 cm (W) x 90 cm (D)

23" (H) x 30" (W) x 35" (D)

Chamber: 38 cm (H) x 43 cm (W) x 83 cm (D) 15" (H) x 17" (W) x 34"(D)

1 unit/case





### Sterilization Equipment

#### Aerators and Abators



#### **Aerator**

The model XL aeration cabinet circulates warm, filtered air in a 135 L (5 cu. ft.) chamber to effectively remove ethylene oxide residual from the goods. Designed for compatibility with  $3M^{\tiny TM}$  Steri-Vac $^{\tiny TM}$  Sterilizers.

- Two selectable aeration temperatures, 37°C (99°F) or 55°C (131°F)
- Modular system for versatility
- Digital display of temperature and elapsed time
- Continuous operation and control via adjustable timer

Cat. No.	Description	Shipping
3M™ Steri-Vac™ XL Aerator		
XL	Aerator	1 unit/case
	Exterior: 62 cm (H) x 78 cm (W) x 93 cm (D)	
	24-3/8" (H) x 30-3/4" (W) x 36-5/8" (D)	
	Chamber: 41 cm (H) x 46 cm (W) x 84 cm (D)	
	16-1/8" (H) x 18-1/8" (W) x 33" (D)	



#### **Ethylene Oxide (EtO) Abator**

The abator removes EtO gas exhausted from  $3M^{TM}$  Steri-Vac<sup>TM</sup> Sterilizers and Aerators, catalytically converting it into  $CO_2$  and water vapor to meet environmental requirements. At normal operating temperatures and concentrations, removal efficiency is  $\geq 99.9\%$ .

- · Works automatically with sterilization system
- · Interfaces with up to two sterilizers
- · Temperature range monitored

Cat. No.	Description	Shipping
3M™ Ethy	lene Oxide (EtO) Abator	
50	Abator	1 unit/case
	80 cm (H) x 82 cm (W) x 105 cm (D)	
	31" (H) x 32" (W) x 41" (D)	

### Gas Cartridges, Monitoring Badges and Printer Papers

#### **Ethylene Oxide Gas Cartridges**

100% EtO is contained in a single-dose cartridge. The cartridge contains 100, 127 or 170 grams (3.5, 4.5 or 6 oz.) of EtO depending on the model Steri-Vac Sterilizer used.

- The system provides an easy, cost-effective and efficient method of gas sterilization
- Precise usage of EtO with single-dose cartridges means no waste of sterilant
- With 3M<sup>™</sup> Steri-Gas<sup>™</sup> Cartridges, operators have no bulky, heavy gas tanks to handle, and no hoses, valves and filters to adjust or replace
- The cartridges do not contain CFCs or HCFCs
- Single dose cannister designed to puncture only when correct vacuum is obtained
- · Custom fill available

Cat. No.	Description	Packaging	Shipping
3M™ Steri-0	Gas <sup>™</sup> Cartridges		
4-100	100 gram Cartridge for Models 4XL and 5XL	12 each/box	4 boxes/case
4-134	127 gram Cartridge for Model 5XL	12 each/box	4 boxes/case
8-170	170 gram Cartridge for Model 8XL	12 each/box	4 boxes/case

#### **Ethylene Oxide Monitoring Badge**

Lightweight personal monitoring badge packed with ethylene oxide-sensitive material. Clips to lapel or pocket near employee breathing zone for passive monitoring of individual during a work shift to meet U.S. mandated requirements. Monitor is returned to 3M support lab with prompt notification of results. Report is mailed to user (or telephoned if permissible exposure limits are exceeded).

Cat. No.	Description	Packaging	Shipping
3M™ Ethyl	ene Oxide (EtO) Monitoring Badge		
3550	EtO Monitoring Badge	_	5 each/case

#### **Printer Papers**

Cat. No.	Description	Shipping
3M™ Steri-\	/ac™ Printer Papers	
1216A	Printer Paper, 60 mm (2-3/8") wide, 5XL (with serial number below 720000)	2 rolls/case
1217	Printer Paper, 79 mm (3-1/8") wide, for 8XL and 5XL (with serial number above 720000)	2 rolls/case







### 3M<sup>™</sup> Attest<sup>™</sup> Rapid Readout Biological Indicators Systems



#### 3M™ Attest™ Rapid Readout Biological Indicator for Steam

The 3M Attest Rapid Readout biological indicator system for use with steam sterilization processing is a proprietary technology that provides the fastest biological indicator reading available. 3M Attest Rapid Readout provides biological indicator readings in as little as three hours post-processing.

- Highly resistant, standardized population of Geobacillus stearothermophilus
- Provides all of the benefits of self-contained technology plus results in three hours for steam
- Reduced incubation can improve inventory levels and increase turns
- · Meets FDA classification as a biological indicator
- Complies with ISO and ANSI/AAMI/ISO 11138:2006, Parts 1 and 3

Cat. No.	Description	Packaging	Shipping
3M™ Attest	™ Rapid Readout Biological Indicator		
1292-S	Steam Biological Indicator (Brown cap)	300 each/box	2 boxes/case



#### 3M™ Attest™ 293 and 290 Auto-readers

3M 293 and 290 Auto-readers are designed to incubate and automatically read the fluorescence of the Attest 1292-S Rapid Readout Biological Indicator in 3 hours. The large capacity 3M Attest 293 Auto-reader is designed specifically for industrial sterilization processing.

- Automatic readout with visual display of indicator result
- The 293 has a 36 well capacity and the 290 has a 12 well capacity
- Developed with significant input from leading medical device manufacturers
- Can help streamline your sterilization monitoring process
- Incubates at  $60^{\circ}\text{C} \pm 2^{\circ}\text{C}$



Cat. No.	Description	Packaging	Shipping
3M™ Attes	t™ Auto-readers		
293	Auto-reader for Steam Large capacity 3M Attest Auto-reader For use with 3M Attest 1292-S Rapid Readout Biological Indicator 36 well capacity	1 unit/box	1 box/case
290	Auto-reader for Steam For use with 3M Attest 1292-S Rapid Readout Biological Indicator 12 well capacity	1 unit/box	1 box/case

## Conventional, Self-contained Biological Indicators for Steam

#### 3M<sup>™</sup> Attest<sup>™</sup> Biological Indicators for Steam Sterilization

3M developed the first self-contained biological indicators to reduce contamination from improper transfer technique.

- Highly resistant, standardized population of Geobacillus stearothermophilus spores
- Conventional self-contained technology benefits:
  - Simplistic culturing
  - Minimizes false positives due to contamination
  - Cost-effective
- Routinely documented for incubation times of less than 7 days (typically 48 hours)
- · Meets FDA classification as a biological indicator
- Complies with ISO and ANSI/AAMI/ISO 11138:2006, Parts 1 and 3

Cat. No.	Description	Packaging	Shipping
3M™ Attest	™ Biological Indicator		
1262-S	Steam Biological Indicator (Brown cap)	300 each/box	2 box/case

### 3M<sup>™</sup> Attest<sup>™</sup> Biological Indicator Incubator

Attest Incubators provide optimal conditions for accurate readout of Attest 1262-S Biological Indicator for Steam

- Convenient 28 vial capacity visible from top
- Dry-block design pre-set to correct incubation temperature

Cat. No.	Description	Packaging	Shipping
Biological	Indicator Incubator		
126	28-vial (2-tier) 56° C, 120 nominal	1 each/box	1 box/case
128	28-vial (2-tier) 56° C (steam) 220/240 volt	1 each/box	1 box/case





### Spore Strips and Media



#### 3M<sup>™</sup> Attest<sup>™</sup> Spore Strips and Culture Media for Steam

Spore strips for industrial sterilization applications

- Highly resistant, standardized population of Geobacillus stearothermophilus spores
- Packaged in peel-open glassine paper pouch
- Comply with ISO and ANSI/AAMI/ISO 11138:2006, Part 1

#### **Culture Media**

Tryptic Soy Broth (TSB) for use with 3M Attest Spore Strips for Steam

- Meets USP recommendation for soybean casein digest broth
- Meets NCCLS sterility and growth promotion guidelines

Cat. No.	Description	Packaging	Shipping
3M™ Attest	Spore Strips for Steam		
1265-S	Biological indicator for Steam	100 each/box	4 boxes/case
Culture Med	dia Test Kit		
1265-M	25 1265-S spore strips and 25 TSB media tubes/box	25 of each/box	1 box/case

### Electronic Test System

(outside of US only and selected countries)

#### 3M™ Electronic Test System for Steam Sterilizers

- Three component monitoring system
  - Sensing unit, data converter, comprehensive data analysis software
- Accurate, reliable, precise and independent data capture
- Fully automatic steam penetration test decision (Pass OR Fail)
- Simple to use and interpret
- Digital, paperless records, distributable over corporate networks

Cat. No.	Description	Packaging	Shipping
3M™ Electi	ronic Test System for Steam Sterilizers		
4008	3M ETS Sensing Unit	1 each/box	1 box/case
4009	3M ETS Data Converter	1 each/box	1 box/case
4010	3M ETS Software	1 each/box	1 box/case



### **Chemical Indicators**





For use in parametric release processes when independent monitoring information is indicated.

- · Indicates sterilization critical parameters have been met
- 1243 Moving front technology provides easy-to-read accept/reject results
- For use with suitable steam sterilization cycles
- Comply with ISO and ANSI/AAMI/ISO 11140:2005, Part 1

Cat. No.	Description	Packaging	Shipping
3M™ Comp	ly™ (SteriGage™) Steam Chemical Integrating Indi	cators	
1243	Steam Chemical Integrator, 500 each/bag Class 5 chemical integrating indicator, moving front for gas mixture sterilizers, 5,1 cm x 1,9	•	2 bags/case
Long Strip	ly <sup>™</sup> Steam Chemical Indicator Strip of chemical indicator ink printed on perforated-st brown. This Class 4 multi-variable indicator react riables.	, , ,	
1250	Steam Chemical Indicator Strip 1,5 cm x 20 cm (5/8" x 8")	240 each/box	1 box/case



#### **Process Indicator Label**

#### **Steam Sterilization Indicator Labels**

Peel-off, pressure-sensitive labels for hard-surfaced items such as unwrapped glass or metal.

Cat. No.	Description	Packaging	Shipping
3M™ Comply™	Chemical Indicator Label		
Color changes	from purple to green during steam sterilization		
00160	Steam — For use in steam sterilizers 1,9 cm x 3,2 cm (3/4" x 1-1/4")	1000 labels/roll 1 roll/box	12 rolls/case
1222CH (Europe only)	Steam Chemical Indicator Label 100 mm x 25 mm process indicator labels applied to clear liner	500 labels/roll 1 roll/box	4 box/case

### **Indicator Tapes**

#### **Indicator Tape for Steam Sterilization**

A reliable steam indicator tape that provides immediate identification of processed goods. The chemical indicator lines will turn dark when exposed to a steam sterilization process.

- Reliable adhesive adheres to packaging material
- 1255 uses high tack adhesive

Cat. No.	Description	Packaging	Shipping
3M™ Comply	™ Indicator Tape for Steam Sterilization*		
1222-6N	1,9 cm x 55 m (3/4" x 60 yds.), Tan	_	24 rolls/case
1255-6N	1,9 cm x 55 m (3/4" x 60 yds.), Blue	_	24 rolls/case







#### **Indicator Tape Dispensers**

Cat. No.	Description	Packaging	Shipping
C22 Heavy	Duty Tape Dispenser		
two rolls of	eavy duty tape dispenser is made of rugge f 3/4-inch wide indicator tape. ted base	ed, high-impact, durable plas	stic. Capable of handling
C22	Heavy duty, holds 2 rolls	_	2 each/case
M920 Defir	nite-length Tape Dispenser		
• Set le	definite-length tape dispenser adjusts to on ngth up to 4 inches (10 cm) one indicator tape roll up to 1 inch (2.5 cm)		1.
M920	Definite length	_	6 each/case



### Monitoring Products for Vapor Hydrogen Peroxide Sterilization

#### **Chemical Indicators**



#### **Vapor Phase Hydrogen Peroxide Chemical Indicator**

An internal pack chemical indicator for use in STERRAD® Sterilization Processes (Models 50, 100S, and 200) or other vapor phase hydrogen peroxide applications when validated by the user. Indicator changes from blue to pink when exposed to vapor hydrogen peroxide. Indicator color does not fade when exposed to light.

Cat. No.	Description	Packaging	Shipping
3M™ Comp	oly™ Gas Plasma Chemical Indicator Strip		
1248	Comply Gas Plasma Chemical Indicator Strip, 2,0 cm x 10 cm (13/16" x 4")	250 each/box	4 boxes/case



#### **Indicator Tape for Vapor Phase Hydrogen Peroxide Processes**

A reliable Vapor Hydrogen Peroxide indicator tape that provides immediate identification of processed goods.

- The chemical indicator will turn from blue to pink when exposed to vapor hydrogen peroxide
- Suitable for use in clean room
- · Reliable adhesive adheres to packaging material and other surfaces
- Backing is cellulose-free

Cat. No.	Description	Packaging	Shipping
3M™ Com	ply™ Gas Plasma Indicator Tape		
1228	1,9 cm x 55 m (3/4" x 60 yds.), White	_	24 rolls/case

### Monitoring Products for Dry Heat Sterilization

### **Biological Indicators**

### 3M<sup>™</sup> Attest<sup>™</sup> Spore Strips and Culture Media for Ethylene Oxide or Dry Heat

Spore strips for industrial sterilization applications

- Highly resistant, standardized population of Bacillus atrophaeus spores
- · Packaged in peel-open glassine paper pouch

#### **Culture Media**

Tryptic Soy Broth (TSB) for use with 3M Attest Spore Strips for Dry Heat or EtO

- Meets USP recommendation for soybean casein digest broth
- Meets NCCLS sterility and growth promotion guidelines

Cat. No.	Description	Packaging	Shipping
3M™ Attest"	* Spore Strips for EtO or Dry Heat		
1263-S	Biological indicator strip for EtO or Dry Heat	100 each/box	4 boxes/case
Culture Med	lia Test Kit		
1263-M	25 1263-S spore strips and 100 TSB media tubes/box	25 of each/box	1 box/case



### Monitoring Products for Dry Heat Sterilization

### **Chemical Indicators**



#### **Dry Heat Chemical Indicator**

Class 4 multi-parameter chemical indicator for dry heat sterilization cycles. Indicator changes color from tan to black within 5 minutes at 170°C (340°F) during exposure to dry heat processing.

• Convenient size

Cat. No.	Description	Packaging	Shipping
3M™ Comp	ly™ Dry Heat Chemical Indicator Strip		
00311	Dry Heat Chemical Indicator Strip 1,5 cm x 11 cm (5/8" x 4-1/4")	250 each/box	10 boxes/case



#### **Sterilization Indicator Labels**

Peel-off, pressure-sensitive labels for hard-surfaced items.

- · Process indicator
- Easy to apply
- The indicator ink on the label turns from tan to black within 5 minutes at 170°C (340°F) to verify exposure during dry heat sterilization.

Cat. No.	Description	Packaging	Shipping
3M™ Compl	ly™ Dry Heat Chemical Indicator Strip		
00170	Dry Heat — For use in dry heat sterilizers 170°C (340°F) or above 1,9 cm x 3,2 cm (3/4" x 1-1/4")	1000 labels/roll 1 roll/box	12 rolls/case

### Other Sterilization Monitoring Products

### Packaging

#### **Sterilization Packages**

Cat. No.	Description	Packaging	Shipping
3M™ Steri-	Lok™ Closure Tape		
1225	White flatback paper tape 2,5 cm x 55 m (1" x 60 yds.)	18 rolls/case	case
3M™ Steri-	Lok™ Peel Open Packages for Ethylene Oxide (E	tO) Sterilization (Film/Tyv	/ek™ Component)
8501	7,9 cm x 20 cm (3-1/8" x 7-7/8")	250 each/box	4 boxes/case
8502	7,9 cm x 24,8 cm (3-7/8" x 9-3/4")	250 each/box	4 boxes/case
8503	14,9 cm x 24,8 cm (5-7/8" x 9-3/4")	250 each/box	4 boxes/case
8504	7,9 cm x 61 cm (3-1/8" x 24")	250 each/box	4 boxes/case
8505	17,8 cm x 33 cm (7" x 13")	250 each/box	2 boxes/case
8506	29,8 cm x 38 cm (11-3/4" x 14-7/8")	250 each/box	2 boxes/case
8507	40 cm x 45,7 cm (15-3/4" x 18")	250 each/box	2 boxes/case
3M™ Steri-	Lok™ Polyethylene Bags for Ethylene Oxide (EtO	) Sterilization	
8250	8,3 cm x 14 cm (3-1/4" x 5-1/2")	100 each/pack	5 packs/case
8251	17,8 cm x 33 cm (7" x 13")	100 each/pack	•
	17,8 cm x 33 cm (7" x 13") 22,9 cm x 43,2 cm (9" x 17")	100 each/pack 100 each/pack	5 packs/case
8252	· · · · · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · · · · ·	5 packs/case 5 packs/case
8252 8253	22,9 cm x 43,2 cm (9" x 17")	100 each/pack	5 packs/case 5 packs/case 5 packs/case
8251 8252 8253 8254 8255	22,9 cm x 43,2 cm (9" x 17") 30,5 cm x 38 cm (12" x 15")	100 each/pack 100 each/pack	5 packs/case 5 packs/case 5 packs/case 5 packs/case 5 packs/case
8252 8253 8254 8255	22,9 cm x 43,2 cm (9" x 17") 30,5 cm x 38 cm (12" x 15") 38 cm x 58,4 cm (15" x 23")	100 each/pack 100 each/pack 100 each/pack 50 each/pack	5 packs/case 5 packs/case 5 packs/case 5 packs/case
8252 8253 8254 8255	22,9 cm x 43,2 cm (9" x 17") 30,5 cm x 38 cm (12" x 15") 38 cm x 58,4 cm (15" x 23") 50,9 cm x 76 cm (20" x 30")	100 each/pack 100 each/pack 100 each/pack 50 each/pack	5 packs/case 5 packs/case 5 packs/case 5 packs/case 5 packs/case
8252 8253 8254 8255 3M <sup>™</sup> Steri-	22,9 cm x 43,2 cm (9" x 17") 30,5 cm x 38 cm (12" x 15") 38 cm x 58,4 cm (15" x 23") 50,9 cm x 76 cm (20" x 30")  Lok™ Polyethylene Tubes for Ethylene Oxide (Etc.)	100 each/pack 100 each/pack 100 each/pack 50 each/pack	5 packs/case 5 packs/case 5 packs/case 5 packs/case 5 packs/case 4 boxes/case
8252 8253 8254 8255 3M <sup>™</sup> Steri- 8203	22,9 cm x 43,2 cm (9" x 17")  30,5 cm x 38 cm (12" x 15")  38 cm x 58,4 cm (15" x 23")  50,9 cm x 76 cm (20" x 30")  Lok** Polyethylene Tubes for Ethylene Oxide (Etc. 7,6 cm x 22,9 m (3" x 75 ft.)	100 each/pack 100 each/pack 100 each/pack 50 each/pack  O) Sterilization 1 each/box	5 packs/case 5 packs/case 5 packs/case 5 packs/case





### Other Sterilization Monitoring Products









	Description	Packaging	Shipping
3M™ Attest	™ Rapid Readout Biological Indicator Test Packs		
	Rapid Readout Test Pack for Steam for steam sterilizers (Geobacillus stearothermophilus) If Attest 1292 Rapid Readout Biological Indicator	25 items/box	25 boxes/case
	Rapid Readout Test Pack for EtO for ethylene oxide (Bacillus atrophaeus) If Attest 1294 Rapid Readout Biological Indicator for EtO	25 items/box	25 boxes/case
3M™ Attest	™ Biological Indicator Test Packs		
,	Test Pack for Steam – Cap color: blue for steam sterilizers (Geobacillus stearothermophilus) A Attest 1262 Steam Indicator	25 items/case	_
1278 Test packs f	Test Pack for EtO – Cap color: green for use with EtO sterilizers (Bacillus atrophaeus)	25 items/case	_
2134MM	Steam Chemical Integrators  Steam chemical integrator  3/4" x 4"  yrators; 250°F/121°C and 270°F/132°C gravity and vacual	250 units/pack num-assisted steam ste	8 packs/case erilization processes
2134MM Steam integ	Steam chemical integrator 3/4" x 4"		
2134MM Steam integ 3M™ Compl 2524MM	Steam chemical integrator 3/4" x 4" rators; 250°F/121°C and 270°F/132°C gravity and vacu	num-assisted steam ste	
2134MM Steam integ 3M** Compl 2524MM EtO integrat	Steam chemical integrator 3/4" x 4"  rators; 250°F/121°C and 270°F/132°C gravity and vacual  ly™ Thermalog™ EtO Chemical Integrators  EtO chemical integrator 3/4" x 4"	num-assisted steam ste	erilization processes
2134MM Steam integ 3M** Compl 2524MM EtO integrat	Steam chemical integrator 3/4" x 4"  rators; 250°F/121°C and 270°F/132°C gravity and vacually Thermalog™ EtO Chemical Integrators  EtO chemical integrator 3/4" x 4"  ors; For use in EtO/HCFC gas mixture sterilization proces	num-assisted steam ste	erilization processes
2134MM Steam integ 3M™ Compl 2524MM EtO integrate 3M™ Compl	Steam chemical integrator 3/4" x 4"  rators; 250°F/121°C and 270°F/132°C gravity and vacually™ Thermalog™ EtO Chemical Integrators  EtO chemical integrator 3/4" x 4"  ors; For use in EtO/HCFC gas mixture sterilization process  ly™ Bowie-Dick™ Type Test Pack  Bowie-Dick Internal Steam	num-assisted steam ste  250 units/pack ses	rilization processes 8 packs/case

### Performance Standards References

#### **Biological Indicators**

International Standards Organization (ISO). ISO 11138 – 1. Sterilization of healthcare products: Biological Indicators – Part 1: General Requirements (2006)

International Standards Organization (ISO). ISO 11138 – 2. Sterilization of healthcare products: Biological Indicators – Part 2: Biological Indicators for ethylene oxide sterilization (2006)

International Standards Organization (ISO). ISO 11138 – 3. Sterilization of healthcare products: Biological Indicators – Part 3: Biological indicators for moist heat sterilization (2006)

International Standards Organization (ISO). ISO 11138 – 4. Sterilization of healthcare products: Biological Indicators – Part 4: Biological indicators for dry heat sterilization (2006)

International Standards Organization (ISO). ISO 14161. Sterilization of healthcare products: Biological Indicators – Guidance for the selection, use and interpretation of results. (2000)

#### **Chemical Indicators**

International Standards Organization (ISO). ISO 11140 – 1. Sterilization of healthcare products: Chemical Indicators – Part 1: General Requirements (2005)

International Standards Organization (ISO). ISO 15882. Sterilization of healthcare products: Chemical Indicators – Guidance for the selection, use and interpretation of results. (2003)

### Index

#### **Index by Catalog Number**

Catalog Number	Product Name	Page Number
00160	3M™ Comply™ Chemical Indicator Label	18
00170	3M™ Comply™ Dry Heat Chemical Indicator Strip	22
0311	3M™ Comply™ Dry Heat Chemical Indicator Strip	22
216A	3M™ Steri-Vac™ Printer Papers	13
217	3M™ Steri-Vac™ Printer Papers	13
222CH	3M™ Comply™ Chemical Indicator Label	18
222-6N	3M™ Comply™ Indicator Tape for Steam Sterilization	19
224-6	3M™ Comply™ Indicator Tape for Ethylene Oxide (Et0) Sterilization	10
225	3M™ Steri-Lok™ Closure Tape	23
228	3M™ Comply™ Gas Plasma Indicator Tape	20
243	3M™ Comply™ (SteriGage™) Steam Chemical Integrating Indicators	18
244A	3M™ Comply™ (SteriGage™) Ethylene Oxide (Et0) Chemical Integrating Indicators	9
248	3M™ Comply™ Gas Plasma Chemical Indicator Strip	20
250	3M™ Comply™ Steam Chemical Indicator Strip	18
251	3M™ Comply™ Ethylene Oxide (Et0) Chemical Indicator Strip	9
255-6N	3M™ Comply™ Indicator Tape for Steam Sterilization	19
26	3M™ Biological Indicator Incubator - 56°C	15
262-S	3M™ Attest™ Bioligical Indicator	15
263-M	Culture Media Test Kit	8
263-M	Culture Media Test Kit	21
263-S	3M™ Attest™ Spore Strips for EtO or Dry Heat	8
263-S	3M™ Attest™ Spore Strips for Et0 or Dry Heat	21
264-S	3M™ Attest™ Biological Indicators for Ethylene Oxide (Et0) Sterilization	7
265-M	Culture Media Test Kit	16
265-S	3M™ Attest™ Spore Strips for Steam	16
27	3M™ Attest™ Biological Indicator Incubator - 37°C	7
276	3M™ Attest™ Biological Indicator Test Packs	24
278	3M™ Attest™ Biological Indicator Test Packs	24
28	3M™ Biological Indicator Incubator - 56°C	15
292-S	3M™ Attest™ Rapid Readout Biological Indicator	14
294-S	3M™ Attest™ Rapid Readout Biological Indicator for Ethylene Oxide	6
296	3M™ Attest™ Rapid Readout Biological Indicator Test Packs	24
290 29	3M <sup>TM</sup> Attest <sup>TM</sup> Biological Indicator Incubator - 37°C	7
298	3M™ Attest™ Rapid Readout Biological Indicator Test Packs	24
134MM	3M™ Comply™ Thermalog™ Steam Chemical Integrators	24
	3M™ Comply™ Thermalog™ Steam Chemical Integrators  3M™ Comply™ Thermalog™ EtO Chemical Integrators	
524MM		24
90	3MTM AttestTM Auto-readers for Steam	14
90G	3MTM AttestTM Auto-readers for Steam	6
93	3MTM AttestTM Auto-readers for Ethylene Oxide	14
93G	3M™ Attest™ Auto-readers for Ethylene Oxide	6
550	3M™ Ethylene Oxide (Et0) Monitoring Badge	13
800	3M ETS Sensing Unit	17
009	3M ETS Data Converter	17
1010	3M ETS Software	17

#### **Index by Catalog Number**

Catalog Number	Product Name	Page Number
4-100	3M™ Steri-Gas™ Cartridges	13
4-134	3M™ Steri-Gas™ Cartridges	13
50	3M™ Ethylene Oxide (EtO) Abator	12
5XL	3M™ Steri-Vac™ 5XL Sterilizer/Aerator	11
8-170	3M™ Steri-Gas™ Cartridges	13
8203	3M™ Steri-Lok™ Polyethylene Tubes for Ethylene Oxide (Et0) Sterilization	23
8205	3M™ Steri-Lok™ Polyethylene Tubes for Ethylene Oxide (Et0) Sterilization	23
8210	3M™ Steri-Lok™ Polyethylene Tubes for Ethylene Oxide (Et0) Sterilization	23
8220	3M™ Steri-Lok™ Polyethylene Tubes for Ethylene Oxide (Et0) Sterilization	23
8250	3M™ Steri-Lok™ Polyethylene Bags for Ethylene Oxide (Et0) Sterilization	23
8251	3M™ Steri-Lok™ Polyethylene Bags for Ethylene Oxide (Et0) Sterilization	23
8252	3M™ Steri-Lok™ Polyethylene Bags for Ethylene Oxide (Et0) Sterilization	23
8253	3M™ Steri-Lok™ Polyethylene Bags for Ethylene Oxide (Et0) Sterilization	23
8254	3M™ Steri-Lok™ Polyethylene Bags for Ethylene Oxide (Et0) Sterilization	23
8255	3M™ Steri-Lok™ Polyethylene Bags for Ethylene Oxide (Et0) Sterilization	23
8501	3M™ Steri-Lok™ Peel Open Packages for Ethylene Oxide (Et0) Sterilization (Film/Tyvek™ Component)	23
8502	3M™ Steri-Lok™ Peel Open Packages for Ethylene Oxide (Et0) Sterilization (Film/Tyvek™ Component)	23
8503	3M™ Steri-Lok™ Peel Open Packages for Ethylene Oxide (Et0) Sterilization (Film/Tyvek™ Component)	23
8504	3M™ Steri-Lok™ Peel Open Packages for Ethylene Oxide (Et0) Sterilization (Film/Tyvek™ Component)	23
8505	3M™ Steri-Lok™ Peel Open Packages for Ethylene Oxide (Et0) Sterilization (Film/Tyvek™ Component)	23
8506	3M™ Steri-Lok™ Peel Open Packages for Ethylene Oxide (Et0) Sterilization (Film/Tyvek™ Component)	23
8507	3M™ Steri-Lok™ Peel Open Packages for Ethylene Oxide (Et0) Sterilization (Film/Tyvek™ Component)	23
8XL	3M™ Steri-Vac™ 8XL Sterilizer/Aerator	11
C22	Heavy Duty Tape Dispenser	10
C22	Heavy Duty Tape Dispenser	19
M920	Definite-length Tape Dispenser	10
M920	Definite-length Tape Dispenser	19
XL	3M™ Steri-Vac™ XL Aerator	12

Call us at Customer Service to place an order: 800-742-1994 (U.S.). Visit our website: **www.3M.com/medicalspecialties** for product and services information, news about conferences we will be attending, new product highlights or to make a direct inquiry.

To have a sales representative contact you, or to request samples, please contact us at 3M HELPLINE 800-228-3957 (U.S.) or for international inquiries, please contact your local country representative. Our 3M Medical Specialties subsidiary contacts are listed below.

#### **International Locations**

Argentina

3M Argentina S.A.C.I.F.I.A. Tel.: 54-11-4339-2400 Fax: 54-11-4339-2640

**Australia** 

3M Health Care Tel: 1 300 363 878 Fax: 1 800 060 888

Belgium

3M Belgium N.V./S.A. Tel: 32-2-722-5111 Fax: 32-2-720-0225

Brazil

3M DO Brasil Ltda. Tel: 55-19-864-7171 Fax: 55-19-864-7848

**Canada** 3M Canada Tel: 800 364-3577 Fax: 800 341-4630

Chile

3M Chile S.A. Tel: (56-2) 4103000 Fax: (56-2) 4103400

China

3M China Ltd. Tel: 86-21-6275-3535 Fax: 86-21-5208-2205

Colombia

3M Colombia S.A. Tel: 57-1 4108585 Fax: 57-1 4161677

**Czech Republic** 

3M Cesko s.r.o. Tel: +420-261 380 111 Fax: +420-261 380 110 Denmark

3M a/s Tel: 45-43-48-0100 Fax: 45-43-96-8596

Dominican Republic 3M Dominicana

Tel: (809) 530-6560 X499 Fax: (809) 530-2960

**Finland** 

Suomen 3M Oy Tel.: 358-9-52-521 Fax: 358-9-512-29-44

**France** 

Laboratoires 3M Santé Tel: 33-1-3031-8376 Fax: 33-1-3031-8378

Germany

3M Medica Tel: 49 2131 14 40 00 Fax: 49 2131 14 49 99

**Hong Kong** 

3M Hong Kong Ltd. Tel: (852) 2806 6111 Fax: (852) 2234 6044

Hungary

3M Hungária Kft. Tel: (36-1) 270-7777 Fax: (36-1) 320-0951

India

3M India Ltd. Tel: 0091-80-2223 1414 Fax: 0091-80-22231450

Indonesia

PT 3M Indonesia Tel: (62-21) 520 3401 Fax: (62-21) 520 3106 Italy

3M Italia SPA Tel: +39 2 7035-2406 Fax: +39 2 7035-2484

Janan

3M Health Care Ltd. Tel: 81-3-3709-8289 Fax: 81-3-3709-8754

Kazakhstan

3M Representation Office Tel: +7 3272 509 944 Fax: +7 3272 509 573

Korea

3M Korea Ltd. Tel: 82-2-3771-4286 Fax: 82-2-786-2825

Malaysia

3M Malaysia Sdn. Bhd. Tel: 6-03-7806 2888 Fax: 6-03-7806 2902

Mexico

3M México, S.A. de C.V. Tel: +52-55-5270-0400 Fax: +52-55-5270-0433

Middle East

3M Gulf Ltd. Tel: 971-4-3670-777 Fax: 971-4-3670-699

The Netherlands

3M Nederland B.V. Tel: 31-715-450-450 Fax: 31-715-450-212

New Zealand 3M New Zealand

Tel: 64-9-444 4760 Fax: 64-9-443 7885 Peru

3M Peru S.A. Tel: 511-224-2728 Fax: 511-224-3171

**Philippines** 

3M Philippines Tel: (632) 8133781 to 95 Fax: (632) 8145873 to 74

**Poland** 

3M Poland Sp.z o.o. Tel: (48 22) 739 60 00 Fax: (48 22) 739 60 01

**Portugal** 

3M Portugal Tel: 351-21-3134500 Fax: 351-21-3134680

**Puerto Rico** 

3M Puerto Rico Tel: 787-620-3000 Fax: 787-620-4612

Russia

3M Russia Tel: 7 (495) 784-74-74

Singapore

3M Technologies (S) Pte. Ltd. Tel: 65-6450-8866 Fax: 65-6455-2130

Fax: 7 (495) 784-74-75

Sri Lanka

3M Lanka (Pvt.) Ltd. Tel: (94) 11 2785701 Fax: (94) 11 4410083-4

Spain

3M España, S.A. Tel: 43-1-321-60-00 Fax: 43-1-321-60-02 Sweden

3M Svenska AB Tel: +46-8-92 21 00 Fax: +46-8-92 24 22

Switzerland & Eastern

Europe

3M (Schweiz) AG Tel: +41 44 724 92 44 Fax: +41 44 724 94 80

Thailand

3M Thailand Ltd. Tel: (66)2 260 8577 Ext.175 Fax: (66)2 261 7535

Taiwan

3M Taiwan Ltd. Tel: 886-2-2704-9011 Fax: 886-2-2706-0355

Turkey

3M Turkey

Tel: 9 (0212) 3507777 Fax: 9 (0212) 2821742

Ukraine

3M Ukraine

Tel: +380 44 490 57 77 Fax: +380 44 490 57 75

United Kingdom & Ireland 3M Health Care Ltd.

Tel: 44-1509-611-611 Fax: 44-1509-237-288

Venezuela

3M Manufacturera Venezuela Tel: 58-241-8391911 Fax: 58-241-8391928

**Viet Nam** 

3M Viet Nam Ltd. Tel: 84 8 640986 Fax: 84 8 641024



**Health Care** 

3M Center, Building 275-4W-02 St. Paul, MN 55144-1000 U.S.A. 1 800 228-3957

www.3m.com/healthcare

3M Canada

Post Office Box 5757 London, Ontario N6A4TI Canada 1 800 563-2921 3M, Attest, Comply, Steri-Gage, Steri-Gas, Steri-Lok and Steri-Vac are trademarks of 3M.

STERRAD is a registered trademark of Johnson & Johnson, Inc. Tyvek is a trademark of Dupont.

Please recycle. Printed in U.S.A.. © 3M 2007. All rights reserved.

