

MX-M Microplate Mixer



Specifications

Speed range	0-1500rpm
Shaking movement	Orbital
Orbit diameter	4.5mm
Max. shaking weight	1lbs
Motor type	Brushless
Motor rating, input	18W
Motor rating, output	10W
Permissible relative humidity	80%
Permissible ambient temperature	5 - 40°C
Voltage	100-240VAC, 50-60Hz
Power	20W
Dimensions	10"L x 6"W x 4"H
Weight	7lbs

Features

- Compact universal small shaker suitable for 96 microtiter plates
- Continuous operation
- Stepless speed regulation of motor
- External rotor brushless motor reduces operating noises and increases motor life
- **2 YEAR WARRANTY**



MX-M with double plate holder accessory



Ordering Information

Cat. No.	Description
82200001	MX-M Microplate Mixer c/w single plate holder

Accessories

18900079	Double microplate holder
18900067	Circular adapter for tube foam inserts
18900020	Foam tube insert, 48 tubes at 6mm dia.
18900021	Foam tube insert, 18 tubes at 10mm dia.
18900022	Foam tube insert, 12 tubes at 12mm dia.
18900023	Foam tube insert, 8 tubes at 16mm dia.
18900024	Foam tube insert, 8 tubes at 20mm dia.