

Specifications	ROCKER MF magnetic filter folder				PALL magnetic filter holder	
						
Model	<i>MF3</i>	<i>MF3a</i>	<i>MF5</i>	<i>MF5a</i>	<b>4242</b>	<b>4238</b>
Stem type	<i>Long</i>	<i>Short</i>	<i>Long</i>	<i>Short</i>	<i>Long</i>	<i>Long</i>
Patent No.	<b>US 9,452,373 B1</b>					
Funnel capacity	300 ml		500 ml		300 ml	500 ml
Connection	Strong magnetic seal				Strong magnetic seal	
Filter diameter	47 mm				47 mm	
Effective filtration area	<b>13.1 cm<sup>2</sup></b>				<b>9.6 cm<sup>2</sup></b>	13.1 cm <sup>2</sup>
Graduation Increment	50 ml				50 ml	
Hose barb	4-8 mm				No	
Adaptive container (short stem limited)	<b>Bottle, Flask ( neck ID &lt; 5 cm)</b>				<b>Flask with side arm</b>	
Autoclavable	YES				YES	
<b>Materials</b>						
Filter funnel	PES (Polyethersulfone)				PPS (Polyphenylsulfone)	
Funnel support base	PES (Polyethersulfone)				PPS (Polyphenylsulfone)	
Membrane support	PES (Polyethersulfone)				PPS (Polyphenylsulfone)	
Lid kit	PP (Polypropylene) construction with PTFE syringe filter				PPS (Polyphenylsulfone) construction with PTFE syringe filter (optional item)	
Hose barb	PES (Polyethersulfone)				No	