

90mm COMPLETE BOTANICAL EXTRACTION FILTER KIT

Sterlitech's 90 mm Complete Botanical Extraction Filter Kit enables easy filtration of batch sizes up to 5 liters. This kit includes everything required to start filtering solvent-extracted compounds of interest from plant-based sources. Perform rough debris removal and bacterial filtration on the same assembly.



FEATURES:

- Clamping mechanism to secure membrane in place
- Removable filter support screen enables easy cleaning and maintenance
- Compatible with polymeric 90 mm disc filters

APPLICATIONS:

- Cold-filtration to remove lipids or waxes
- Clarification or plant debris removal
- Bacterial removal (using 0.2 or 0.45 μm PES filters)
- Batch processing with manifolds
- · Ethanol washing

ORDERING INFO:

SKU: 2050001

Shipping Dimensions:

Box One: 25 x 21 x 21 in. - 36 lbs

COMPONENTS:

- 5L side arm flask
- Glass filter holder with PTFE-coated stainless steel support screen
- PILOT 5000 vacuum pump
 - 110/220V
 - UL/CE listed
- EPDM rubber tubing
- Membrane packs:
 - Coarse cellulose filter paper (19-26 μm)
 - Medium cellulose filter paper (6 μm)
 - Fine membrane filters (0.45 µm PES)